

NEW

Smart-Caller
Keeping You Informed

cura1

DOOR ALERT SYSTEM

Helping to care.



Door Alert immediately alerts caregivers when a monitored person exits a protected door or exit. Alternatively Door Alert can silently lock a door, prevent an automatic door from opening or disable a lift for a predetermined time, thereby reducing the need for intervention by caregivers.

SMART-CALLER P/L

03 9588 0833 sales@smartcaller.com.au



This product is approved
by the Therapeutic
Goods Administration

OUR GUARANTEE: All products supplied with full warranty and guaranteed to be free from defects in material and workmanship.

Door Alert

Door Alert can be configured to work in a variety of ways:

- To activate an alarm when a monitored person approaches a protected exit.
- To alert staff (via an optional nurse call connection) when a monitored person approaches a protected exit.
- To secure an exit when a monitored person approaches.



wrist/ankle transmitter


reset transmitter



door monitor



system tester (optional)



- ▷ Low cost
- ▷ Unobtrusive design
- ▷ Fully supervised wireless operation
- ▷ Low cost waterproof transmitters
- ▷ Easy installation and range setting
- ▷ Monitor up to 1000 residents and up to 16 exit points with Care Station
- ▷ Individually identifies each resident
- ▷ Resident transmitters can be switched on/off by authorised persons to conserve battery life
- ▷ Inbuilt audible and visual alerts
- ▷ Alarm output to trigger external alarm devices or nurse call system
- ▷ Low battery indication for all transmitters
- ▷ Magnetic door contact input to inhibit alarm activation if the door is closed
- ▷ Can control 3rd party equipment i.e. door locks, elevators, sliding doors etc.

Visit our web site for more solutions and comprehensive Catalogue

www.smartcaller.com.au